# INSTRUCTIONS FOR USE

Medical compression stockings

- |   |   |   |
|---|---|---|
| <b>(EN)</b> Medical compression stockings     | <b>(NL)</b> Medische compressiekousen                 | <b>(SL)</b> Medicinske kompresijske nogavice  |
| <b>(DE)</b> Medizinische Kompressionsstrümpfe | <b>(PT)</b> Meias de compressão medicinais            | <b>(RU)</b> Медицинские компрессионные чулки  |
| <b>(FR)</b> Bas médicaux de compression       | <b>(PL)</b> Medyczne pończosznicze wyroby kompresyjne | <b>(IS)</b> Teygjusokkar til lækninga         |
| <b>(ES)</b> Medias de compresión médica       | <b>(HU)</b> Kompressziós térdharisnya                 | <b>(AR)</b> جوارب طبية ضاغطة                  |
| <b>(SV)</b> Medicinska kompressionsstrumpor   | <b>(IT)</b> Calze compressive medicali                | <b>(TR)</b> Tibbi kompresyon çorapları        |
| <b>(NO)</b> Medisinske kompresjonsstrømper    | <b>(FI)</b> Lääkinnälliset kompressiosukat            | <b>(HE)</b> גרבי לחץ רפואיים                  |
| <b>(DA)</b> Medicinsk kompressionsstrømpe     | <b>(RO)</b> Ciorapi compresivi medicinali             | <b>(CZ)</b> Zdravotnické kompresivní punčochy |

(EN) FIBER CONTENT

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(FR) COMPOSITION

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(NO) FIBERINNHOLD

(DA) FIBERINDHOLD

(NL) VEZELGEHALTE

(PT) TEOR DE FIBRA

(PL) SKŁAD

(HU) ROSTTARTALOM

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(RO) CONTIN FIBRE

(SL) VSEBNOST VLAKEN

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## Dear Customer,

Thank you for choosing JOBST®, the medical compression stockings that incorporates gradient compression for the management of venous and lymphatic diseases.

JOBST® offers you a wide selection of products for these indications – either as made-to-measure or ready-to-wear compression stockings.

### 1. Application

#### Tips for Wearing JOBST® Medical LegWear:

To help your legs look and feel good, wear your JOBST® Medical LegWear regularly and consistently. Ideally you should put on your compression stocking in the morning, immediately after getting up. You may experience some difficulty putting on the stocking in the beginning but this will become easier after a number of days.

When applying the compression stockings, please do not pull on the edge of the stockings instead follow these recommendations:

Turn the compression stocking inside out all the way to the heel. Pull apart the folded-in foot part using both hands and pull over the foot all the way to the heel.

Gradually pull the compression stocking over instep and heel.

Gradually pull the compression stocking upward by returning the stocking to a right side out position piece-by-piece.

Distribute the stocking evenly by massaging with both hands until the stocking fits and does not form creases. It is particularly helpful to wear rubber gloves for this step.

For a Waist High garment continue working the garment up to the middle of your thigh. Stand up. With both hands, grab a handful of fabric and pull it up to your waist.

Adjust your support so that it is snug at your midsection. This will keep it from sliding down. Pull the midsection to its proper height.

Slip your hands, palm out, into the back and lift the fabric up and away from your skin. Let it fall back to fit snugly over the buttocks.

### Special notes:

- Caution: long fingernails and jewelry may damage the stockings when putting them on. For easier application use rubber gloves or a donning aid.
- Please use detailed measuring tables from our JOBST® catalogs.
- Do not bring the compression stockings into contact with greasy ointments, body lotions or similar substances unless recommended by JOBST®.
- Loose threads or loops should not be cut off, since this may cause holes or runs to form in the stockings.
- The product should be used by a single patient.

### 2. Care instructions

JOBST® compression stockings should, if at all possible, be washed daily, or at least every other day.

For this purpose, you may use a regular mild detergent. Do not use laundry additives such as fabric softener, optical brightener or stain remover. Rinse your compression stockings well.

Roll up the wet compression stockings in a towel after washing and press out excess moisture. Do not wring. Do not dry the compression stockings on a radiator or in direct sunlight. For tumble dryer compatibility please check your compression stocking's product label.

### Silicone band care:

When needed, hand wash silicone band to remove excess creams, lotions, hair and other foreign material.

### 3. Indications

Every JOBST® Medical LegWear product has a specific compression level designed to help specific venous and lymphatic conditions. The compression level is printed on the package label.

The required compression class must be determined by your physician.

Medical compression garments are used for the management of the following indications:

#### Compression Class 1 (CCL 1)

- Tired, heavy, aching legs
- Minor varices without oedema
- Mild swelling of the feet, ankles and legs
- Onset of pregnancy-related varices
- Controlled cardiac oedema
- Arthritic patients

- Primary or secondary lymphoedema in International Society of Lymphology (ISL) stages 0-II with no or minimal shape distortion (or lymphoedema ISL stage III when higher compression is not tolerated)

## Compression Class 2 (CCL 2)

- Tired, heavy, aching legs
- Moderate to severe varices with and without oedema
- Moderate to severe pregnancy-related varices with and without oedema
- Moderate oedema
- Chronic venous insufficiency (CVI) C3-C5 according to CEAP: venous oedema, skin alterations and healed ulcus cruris
- Helps to prevent the recurrence of venous ulcers
- For adjunct management of superficial thrombophlebitis/SVT
- Post venous surgery
- Post sclerotherapy
- Post-traumatic and post-operative oedema
- Arthritic patients with reasonable dexterity
- Primary or secondary lymphoedema in ISL stages II-III with some shape distortion
- Phlebolympoedema
- Lipoedema

## Compression Class 3 (CCL 3)

- Severe varices with or without oedema
- Severe oedema
- Chronic venous insufficiency (CVI) C3-C6 according to CEAP: venous oedema, skin alterations, healed ulcus cruris and ulcus cruris
- To prevent recurrence of venous ulcers
- Post-thrombotic syndrome (PTS)
- For adjunct management of superficial thrombophlebitis / SVT
- Post venous surgery
- Post sclerotherapy
- Active patients with lymphoedema in ISL stage III
- Risk of rapid oedema rebound
- Severe lipoedema

## Compression Class 4 (CCL 4)

- Severe post-thrombotic syndrome (PTS)
- Severe lymphoedema in ISL stage IV
- Severe oedema of non-arterial origin

## 4. Contraindications

JOBST® Medical LegWear should not be worn in cases of:

- Advanced arterial insufficiency including ischemia
- Uncontrolled congestive heart failure
- Untreated septic phlebitis
- Phlegmasia coerulea dolens
- Incompatibility to fabric

## 5. Cautions

- Weeping dermatosis
- Cutaneous infections
- Open wounds must be covered with an

- appropriate dressing before the compression garment is put on.
- Severely compromised skin sensibility and impaired sensitivity of the limb
- Advanced peripheral neuropathy
- Immobility (confinement to bed). Not recommended to wear during bed rest.
- Conditions in which increased venous and lymphatic return is not desired.
- Appropriate hygiene and skin care should be applied.
- Appropriate footwear must be worn when walking to prevent falls.
- Ill-fitting compression garments may cause constriction leading to e.g. skin necrosis, pressure ulcers and pressure related peripheral nerve damage.
- The attending physician decides whether compression in general can be applied to patients suffering from the following diseases: rheumatoid arthritis, complex regional pain syndrome, CRPS (M. Sudeck), malignant lymphedema, gangrene.
- Medical compression garments should only be used after consultation with your doctor or therapist and checked regularly by them on the appropriateness of your therapy, also considering age and health status.
- Patient measurements should be taken by a trained person to determine the correct product size needed and to decide whether a made-to-measure or ready-to-wear garment is required.
- When re-ordering a garment after 6 months we recommend remeasuring the patient as the limb circumference might have changed.
- If you notice discomfort, pain, numbness, changes in skin color, appearance of new wounds or swelling during wear, stop using the product and contact your healthcare professional.
- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of your state
- Caution applicable for made-to-measure Products: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

## 6. Directions for storage / Guarantee

A maximum storage period of up to 36 months applies for serial products. This is followed by a wearing period of not more than 6 months for the compression stocking. The expiry date for storage is printed on the package label with the symbol of an hourglass. Made-to-measure products are products intended for immediate use and no expiry date for storage is indicated on these. Due to wear or slackening of the material, the medical efficacy can be guaranteed for a wearing period of six months at most both for serially produced sizes and made-to-measure products. This assumes correct handling (e.g. product care, donning and doffing).

Please protect your compression garments from direct exposure to sunlight, heat, and moisture. We recommend storage at room temperature.

If the vascular disease persists after 6 months, please consult your attending physician concerning the possible need for a new compression stocking.

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
Thank you for your trust.

Discover more about healthy legs and compression products on:

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Australian Distributor: BSN medical (Aust.) Pty Ltd

L6, 211 Wellington Road, Mulgrave 3170 Victoria Australia

		Artwork ID: JBR14543.01	Review Loop:  3	Printing Colours:   Pantone 281 C
Project Title: C2365_702 JOBST IFU Lower Extremities EMM				
Contact: Artwork.Management@essity.com				
Contact Address: BSN medical GmbH Quickbornstrasse24 20253 Hamburg   Germany				
Ref Number(s): 34802-99432-00				
Languages: GBR, DEU, FRA, ESP, SWE, NOR, DNK, NLD, PRT, POL, HUN, ITA, FIN, ROU, SVN, RUS, ISL, SAU, TUR, ISR, CZE			Comment on Diecut: NoteLOT and expiry-date symbol will be part of variable coding.	
			Date of Review Loop: 07.10.2019	
			Replaces Artwork ID: -JBR12070.01	
			Technical Drawing for print ID: -	

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BSN-JOBST GmbH